



COQ7 Antibody

Product Code	CSB-PA005806GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99807
Immunogen	Human COQ7
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	coenzyme Q7 homolog, ubiquinone (yeast);COQ7;CAT5;CLK-1;CLK1 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	COQ7
Target Details	This protein is similar to a mitochondrial di-iron containing hydroxylase in Saccharomyces cerevisiae that is involved with ubiquinone biosynthesis. Mutations in the yeast gene lead to slower development and longer life span. Transcript variants have been described for this gene, but their full-length sequences have not been determined.
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.