



DHRS9 Antibody

Product Code	CSB-PA006867GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BPW9
Immunogen	Human DHRS9
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	dehydrogenase/reductase (SDR family) member 9;DHRS9;3ALPHA-HSD;RDH15;RDHL;RDHTBE;RETSR8;SDR9C4 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	DHRS9
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.