



DNMT3A Antibody

Product Code	CSB-PA007083GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y6K1
Immunogen	Human DNMT3A
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB; Recommended dilution: WB:1:1000-1:4000
Form	Liquid
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	DNA (cytosine-5-)-methyltransferase 3 alpha;DNMT3A;DNMT3A2;M.HsaIIIA ;
Immunogen Species	Homo sapiens (Human)
Target Names	DNMT3A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.