



# DTD1 Antibody

<b>Product Code</b>	CSB-PA007214GA01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q8TEA8
<b>Immunogen</b>	Human DTD1
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human,Mouse,Rat
<b>Tested Applications</b>	ELISA,WB
<b>Storage Buffer</b>	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
<b>Purification Method</b>	Antigen Affinity purified
<b>Isotype</b>	IgG
<b>Alias</b>	D-tyrosyl-tRNA deacylase 1 homolog (S. cerevisiae);DTD1;C20orf88;DUEB;HARS2;MGC119131;MGC41905;bA379J5.3; bA555E18.1;pqn-68 ;
<b>Product Type</b>	Purified Rabbit Anti human PolyClonal Antibody
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	DTD1
<b>Target Details</b>	This protein is similar in sequence to histidyl-tRNA synthetase, which hydrolyzes D-tyrosyl-tRNA(Tyr) into D-tyrosine and free tRNA(Tyr). The encoded protein is found in the cytoplasm.
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.