



DUS3L Antibody

Product Code	CSB-PA007236GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96G46
Immunogen	Human DUS3L
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	dihydrouridine synthase 3-like (S. cerevisiae);DUS3L;DUS3;FLJ13896 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	DUS3L
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.