



EGFL7 Antibody

Product Code	CSB-PA007476GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UHF1
Immunogen	Human EGFL7
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	EGF-like-domain, multiple 7;EGFL7;MGC111117;RP11-251M1.2;VE-STATIN;ZNEU1 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	EGFL7
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.