



EIF2A Antibody

Product Code	CSB-PA007509GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BY44
Immunogen	Human EIF2A
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	ELISA, WB, IHC, IF
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	eukaryotic translation initiation factor 2A, 65kDa; EIF2A; CDA02; EIF-2A; MST089; MSTP004; MSTP089 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	EIF2A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.