



# EPB41L2 Antibody

<b>Product Code</b>	CSB-PA007708GA01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	O43491
<b>Immunogen</b>	Human EPB41L2
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human,Mouse,Rat
<b>Tested Applications</b>	ELISA,WB
<b>Storage Buffer</b>	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
<b>Purification Method</b>	Antigen Affinity purified
<b>Isotype</b>	IgG
<b>Alias</b>	erythrocyte membrane protein band 4.1-like 2;EPB41L2;4.1-G;DKFZp781D1972;DKFZp781H1755 ;
<b>Product Type</b>	Purified Rabbit Anti human PolyClonal Antibody
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	EPB41L2
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.