



NHLRC1 Antibody

Product Code	CSB-PA015790GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6VVB1
Immunogen	Human NHLRC1
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	NHL repeat containing 1;NHLRC1;EPM2A;EPM2B;MALIN;MGC119262;MGC119264;MGC119265;bA2 04B7.2 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	NHLRC1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.