



PCNXL2 Antibody

Product Code	CSB-PA017629GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A6NKB5
Immunogen	Human PCNXL2
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	pecanex-like 2 (Drosophila);PCNXL2;FLJ11383;KIAA0435 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PCNX2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.