



PDS5A Antibody

Product Code	CSB-PA017742GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q29RF7
Immunogen	Human PDS5A
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IF
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	PDS5, regulator of cohesion maintenance, homolog A (S. cerevisiae);PDS5A;DKFZp686B19246;FLJ41012;KIAA0648;MGC131948;MGC161503;PIG54;SCC-112 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PDS5A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.