



POLR1A Antibody

Product Code	CSB-PA018322GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95602
Immunogen	Human POLR1A
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	polymerase (RNA) I polypeptide A, 194kDa;POLR1A;DKFZp586M0122;FLJ21915;MGC87965;RPA1;RPO1-4 ;RPA194
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	POLR1A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.