



RPAP1 Antibody

Product Code	CSB-PA020094GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BWH6
Immunogen	Human RPAP1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	RNA polymerase II associated protein 1;RPAP1;DKFZp727M111;FLJ12732;KIAA1403;MGC858 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	RPAP1
Target Details	The function of this gene is unknown.
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.