



RSAD2 Antibody

Product Code	CSB-PA020536GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WXG1
Immunogen	Human RSAD2
Species Reactivity	Human, Mouse, Rat
Tested Applications	ELISA, IF
Form	Liquid
Storage Buffer	PBS with 0.02% sodium azide and 50% glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	RSAD2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.