



SDCCAG8 Antibody

Product Code	CSB-PA020899GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86SQ7
Immunogen	Human SDCCAG8
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat,Zebrafish
Tested Applications	ELISA,WB,IF
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	serologically defined colon cancer antigen 8;SDCCAG8;CCCAP;HSPC085;NY-CO-8 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SDCCAG8
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.