



SH3PXD2A Antibody

Product Code	CSB-PA021239GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5TCZ1
Immunogen	Human SH3PXD2A
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	SH3 and PX domains 2A;SH3PXD2A;FISH;SH3MD1 ;
Product Type	Purified Rabbit Anti Human,Mouse PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SH3PXD2A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.