



SHCBP1 Antibody

Product Code	CSB-PA021258GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NEM2
Immunogen	Human SHCBP1
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	SHC SH2-domain binding protein 1; SHCBP1; FLJ22009; MGC26900; PAL ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SHCBP1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.