



SLC10A3 Antibody

Product Code	CSB-PA021375GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P09131
Immunogen	Human SLC10A3
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	solute carrier family 10 (sodium/bile acid cotransporter family), member 3; SLC10A3; DXS253E; P3 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SLC10A3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.