



SLC17A5 Antibody

Product Code	CSB-PA021421GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NRA2
Immunogen	Human SLC17A5
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	solute carrier family 17 (anion/sugar transporter), member 5;SLC17A5;AST;FLJ22227;FLJ23268;ISSD;NSD;SD;SIALIN;SIASD;SLD ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SLC17A5
Target Details	This gene encodes a membrane transporter that exports free sialic acids that have been cleaved off of cell surface lipids and proteins from lysosomes. Mutations in this gene cause sialic acid storage diseases, including infantile sialic acid storage disorder and and Salla disease, an adult form.
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.