



SLC4A8 Antibody

Product Code	CSB-PA021673GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q2Y0W8
Immunogen	Human SLC4A8
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	solute carrier family 4, sodium bicarbonate cotransporter, member 8;SLC4A8;DKFZp761B2318;FLJ46462;NBC3 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SLC4A8
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.