



SMARCB1 Antibody

Product Code	CSB-PA021805GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q12824
Immunogen	Human SMARCB1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	SWI/SNF related, matrix associated, actin dependent regulator of chromatin, subfamily b, member 1;SMARCB1;BAF47;INI1;RDT;SNF5;SNF5L1;Sfh1p;Snr1;hSNFS ;
Product Type	Purified Rabbit Anti Human,Mouse PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SMARCB1
Target Details	This protein is part of a complex that relieves repressive chromatin structures, allowing the transcriptional machinery to access its targets more effectively. The encoded nuclear protein may also bind to and enhance the DNA joining activity of HIV-1 integrase. This gene has been found to be a tumor suppressor, and mutations in it have been associated with malignant rhabdoid tumors. Two transcript variants encoding different isoforms have been found for this gene.
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.