



SP6 Antibody

Product Code	CSB-PA022447GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q3SY56
Immunogen	Human SP6
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	Sp6 transcription factor;SP6;EPFN;EPIPROFIN;KLF14;MGC119662;MGC119663;MGC119664;MGC119665 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SP6
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.