



SPAG5 Antibody

Product Code	CSB-PA022466GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96R06
Immunogen	Human SPAG5
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA,WB,IF,IP
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	sperm associated antigen 5;SPAG5;DEEPEST;MAP126;hMAP126 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SPAG5
Target Details	This gene encodes a protein associated with the mitotic spindle apparatus. The encoded protein may be involved in the functional and dynamic regulation of mitotic spindles.
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.