



SUZ12 Antibody

Product Code	CSB-PA022976GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15022
Immunogen	Human SUZ12
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	suppressor of zeste 12 homolog (Drosophila);SUZ12;CHET9;JJAZ1;KIAA0160 ;
Product Type	Purified Rabbit Anti Human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SUZ12
Target Details	This zinc finger gene has been identified at the breakpoints of a recurrent chromosomal translocation reported in endometrial stromal sarcoma. Recombination of these breakpoints results in the fusion of this gene and JAZF1. This protein contains a zinc finger domain in the C terminus of the coding region.
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.