



TLR5 Antibody

Product Code	CSB-PA023604GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60602
Immunogen	Human TLR5
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	toll-like receptor 5; TLR5; FLJ10052; MGC126430; MGC126431; SLEB1; TIL3 ;
Product Type	Purified Rabbit Anti Homo sapiens PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TLR5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.