



TRPC4 Antibody

| | |
|----------------------------|-------------------------------------------------------------------------------------------------------------------------------|
| Product Code | CSB-PA025078GA01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q9UBN4 |
| Immunogen | Human TRPC4 |
| Raised In | Rabbit |
| Species Reactivity | Human,Mouse,Rat |
| Tested Applications | ELISA,WB |
| Storage Buffer | PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles. |
| Purification Method | Antigen Affinity purified |
| Isotype | IgG |
| Alias | transient receptor potential cation channel, subfamily C, member 4;TRPC4;HTRP4;MGC119570;MGC119571;MGC119572;MGC119573;TRP4 ; |
| Product Type | Purified Rabbit Anti human PolyClonal Antibody |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | TRPC4 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |