



TRPV6 Antibody

Product Code	CSB-PA025098GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H1D0
Immunogen	Human TRPV6
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IF
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	transient receptor potential cation channel, subfamily V, member 6;TRPV6;ABP/ZF;CAT1;CATL;ECAC2;HSA277909;LP6728;ZFAB ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TRPV6
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.