



TSPAN7 Antibody

Product Code	CSB-PA025165GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P41732
Immunogen	Human TSPAN7
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB
Form	Liquid
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TSPAN7
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.