



UBP1 Antibody

Product Code	CSB-PA025508GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NZI7
Immunogen	Human UBP1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	upstream binding protein 1 (LBP-1a);UBP1;DKFZp686L1745;LBP-1B;LBP-1a;LBP1A;LBP1B ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	UBP1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.