



UGT2B7 Antibody

Product Code	CSB-PA025598GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P16662
Immunogen	Human UGT2B7
Raised In	Rabbit
Species Reactivity	Human,Rat
Tested Applications	ELISA,WB,IHC,IF
Form	Liquid
Storage Buffer	PBS with 0.02% sodium azide and 50% glycerol, pH 7.3.
Purification Method	Antigen Affinity purified
Isotype	IgG
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	UGT2B7
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.