



ZCCHC11 Antibody

Product Code	CSB-PA026374GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5TAX3
Immunogen	Human ZCCHC11
Raised In	Rabbit
Species Reactivity	Human,Mouse
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	zinc finger, CCHC domain containing 11;ZCCHC11;PAPD3 ;
Product Type	Purified Rabbit Anti Human,Mouse PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ZCCHC11
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.