



# Recombinant Human Cytochrome P450 2C19 (CYP2C19)

<b>Product Code</b>	CSB-YP006417HU
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P33261
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	RGKLPPGPTPLPVIGNILQIDIKDVSKSLTNLSKIYGPVFTLYFGLERMVVLHGYYE VVKEALIDLGEEFSGRGRHFPLAERANRGGFIVFSNGKRWKEIRRFSLMTRLRNF GMGKRSIEDRVQEEARCLVEELRKTASPCDPTFILGCAPCNVICSIIIFQKRFD YKDQQFLNLMEKLNENIRIVSTPWIQICNNFPTIIDYFPGTHNKLLKNLAFMESDI LEKVKEHQESMDINNPRDFIDCFLIKMEKEKQNNQSEFTIENLVITAADLLGAGT ETTSTTLRYALLLLLKHPEVTAKVQEEIERVVGRNRSPCMQDRGHMPYTDVV HEVQRYIDLIPISLPHAVTCDVKFRNYLIPKGTTLTSLTSVLHDNKEFPNPEMF DPRHFLDEGGNFKKSNYFMPFSAGKRICVGEGLARMELFLFTFILQNFNLKSL IDPKDLDTTPVVNGFASVPPFYQLCFIPV
<b>Source</b>	Yeast
<b>Target Names</b>	CYP2C19
<b>Protein Names</b>	Recommended name: Cytochrome P450 2C19 EC= 1.14.13.-Alternative name(s): (R)-limonene 6-monooxygenase EC= 1.14.13.80 (S)-limonene 6- monooxygenase EC= 1.14.13.48 (S)-limonene 7-monooxygenase EC= 1.14.13.49 C
<b>Expression Region</b>	26-490
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Target Details</b>	This gene encodes a member of the cytochrome P450 superfamily of enzymes. The cytochrome P450 proteins are monooxygenases which catalyze many reactions involved in drug metabolism and synthesis of cholesterol, steroids and other lipids. This protein localizes to the endoplasmic reticulum and is known to metabolize many xenobiotics, including the anticonvulsive drug mephenytoin, omeprazole, diazepam and some barbiturates. Polymorphism within this gene is associated with variable ability to metabolize mephenytoin, known as the poor metabolizer and extensive metabolizer phenotypes. The gene is located within a cluster of cytochrome P450 genes on chromosome 10q24.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.