



Cortisol Monoclonal Antibody

Product Code	CSB-MA00371110m
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Immunogen	Cortisol-BSA conjugate
Raised In	mouse
Specificity	No significant cross-reactivity or interference was observed
Tested Applications	ELISA
Form	liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, protein G purified
Isotype	IgG1
Clonality	monoclonal
Alias	Hydrocortisone
Product Type	Small Molecules Antibodies
Target Names	COR
Accession NO.	9E1C1