




CDK4 Antibody

Product Code	CSB-PA11659A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P11802
Immunogen	Recombinant Human Cyclin-dependent kinase 4 protein (2-303AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cyclin-dependent kinase 4 (EC 2.7.11.22) (Cell division protein kinase 4) (PSK-J3), CDK4
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	CDK4

Image



Immunohistochemistry of paraffin-embedded EC109 cells using CSB-PA11659A0Rb at dilution of 1:50

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.