




# PSMB1 Antibody

<b>Product Code</b>	CSB-PA01395A0Rb
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P20618
<b>Immunogen</b>	Recombinant Human Proteasome subunit beta type-1 protein (29-241AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Proteasome subunit beta type-1 (EC 3.4.25.1) (Macropain subunit C5) (Multicatalytic endopeptidase complex subunit C5) (Proteasome component C5) (Proteasome gamma chain), PSMB1, PSC5
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Cell Biology
<b>Target Names</b>	PSMB1

## Image



### Western blot

All lanes: Proteasome subunit beta 1 antibody at 2µg/ml  
 Lane 1: EC109 whole cell lysate  
 Lane 2: 293T whole cell lysate  
 Secondary  
 Goat polyclonal to rabbit IgG at 1/15000 dilution  
 Predicted band size: 27 kDa  
 Observed band size: 27 kDa

## Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.