




PSMB1 Antibody

Product Code	CSB-PA01395A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P20618
Immunogen	Recombinant Human Proteasome subunit beta type-1 protein (29-241AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Proteasome subunit beta type-1 (EC 3.4.25.1) (Macropain subunit C5) (Multicatalytic endopeptidase complex subunit C5) (Proteasome component C5) (Proteasome gamma chain), PSMB1, PSC5
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	PSMB1

Image



Western blot

All lanes: Proteasome subunit beta 1 antibody at 2μg/ml

Lane 1: EC109 whole cell lysate

Lane 2: 293T whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/15000 dilution

Predicted band size: 27 kDa

Observed band size: 27 kDa