



RPL17 Antibody, Biotin conjugated

Product Code	CSB-PA02505D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P18621
Immunogen	Recombinant Human 60S ribosomal protein L17 protein (2-184AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	60S ribosomal protein L17 (60S ribosomal protein L23) (Large ribosomal subunit protein uL22) (PD-1), RPL17
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	RPL17
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.