



CTSK Antibody

Product Code	CSB-PA15119A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P43235
Immunogen	Recombinant Human Cathepsin K protein (115-329AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cathepsin K (EC 3.4.22.38) (Cathepsin O) (Cathepsin O2) (Cathepsin X), CTSK, CTSO CTSO2
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	CTSK
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.