



Abcb4 Antibody

Product Code	CSB-PA001050LA01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P21440
Immunogen	Recombinant Mouse Phosphatidylcholine translocator ABCB4 protein (352-708AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Phosphatidylcholine translocator ABCB4 (ATP-binding cassette sub-family B member 4) (Multidrug resistance protein 2) (Multidrug resistance protein 3) (EC 3.6.3.44) (P-glycoprotein 2) (P-glycoprotein 3), Abcb4, Mdr2 Pgy-2 Pgy2
Immunogen Species	Mus musculus (Mouse)
Research Area	Others
Target Names	Abcb4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.