





MDK Antibody

Product Code	CSB-PA10574A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P21741
Immunogen	Recombinant Human Midkine protein (21-143AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Midkine (MK) (Amphiregulin-associated protein) (ARAP) (Midgestation and kidney protein) (Neurite outgrowth-promoting factor 2) (Neurite outgrowth-promoting protein), MDK, MK1 NEGF2
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	MDK

Image



Immunofluorescent analysis of HepG2 cells using CSB-PA10574A0Rb at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)



Immunohistochemistry of paraffin-embedded human small intestine tissue using CSB-PA10574A0Rb at dilution of 1:100

**Usage**

For Research Use Only. Not for use in diagnostic or therapeutic procedures.