



FABP3 Antibody, Biotin conjugated

Product Code	CSB-PA07764D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P05413
Immunogen	Recombinant Human Fatty acid-binding protein, heart protein (2-133AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Fatty acid-binding protein, heart (Fatty acid-binding protein 3) (Heart-type fatty acid-binding protein) (H-FABP) (Mammary-derived growth inhibitor) (MDGI) (Muscle fatty acid-binding protein) (M-FABP), FABP3, FABP11 MDGI
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	FABP3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.