



Cholic acid Monoclonal Antibody

Product Code	CSB-MA005511I0m
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Immunogen	Cholic acid-BSA conjugate
Raised In	mouse
Tested Applications	ELISA
Form	liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%?protein G purified
Isotype	IgG1
Clonality	monoclonal
Alias	Cholate, cholanic acid,cholalin, colalin, cholalic acid, 3alpha,7alpha,12alpha-trihydroxy-5beta-cholanicacid, cholsaeure
Product Type	Small Molecules Antibodies
Target Names	CA
Accession NO.	1C1F1