



Cholic acid Monoclonal Antibody

| | |
|----------------------------|--|
| Product Code | CSB-MA005511I0m |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Immunogen | Cholic acid-BSA conjugate |
| Raised In | mouse |
| Tested Applications | ELISA |
| Form | liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95% protein G purified |
| Isotype | IgG1 |
| Clonality | monoclonal |
| Alias | Cholate, cholanic acid, cholalin, colalin, cholalic acid, 3alpha,7alpha,12alpha-trihydroxy-5beta-cholanic acid, cholsaeure |
| Product Type | Small Molecules Antibodies |
| Target Names | CA |
| Accession NO. | 1C1F1 |