



AMA-1 Antibody

| | |
|----------------------------|--|
| Product Code | CSB-PA338457HA01PLT |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P22621 |
| Immunogen | Recombinant Plasmodium falciparum Apical membrane antigen 1 protein (25-546AA) |
| Raised In | Rabbit |
| Species Reactivity | Plasmodium falciparum |
| Tested Applications | ELISA;Not yet tested in other applications. |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | 0.1mol/L NaCl, 15mmol/L NaN ₃ , PH7.2 |
| Purification Method | Caprylic Acid Ammonium Sulfate Precipitation purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Merozoite surface antigen,AMA-1 |
| Immunogen Species | Plasmodium falciparum (isolate FC27 / Papua New Guinea) |
| Target Names | AMA-1 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |