



CASR Antibody, HRP conjugated

Product Code	CSB-PA004558EB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P41180
Immunogen	Recombinant Human Extracellular calcium-sensing receptor protein (20-612AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Extracellular calcium-sensing receptor (CaR) (CaSR) (hCasR) (Parathyroid cell calcium-sensing receptor 1) (PCaR1), CASR, GPRC2A PCAR1
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	CASR