



CASR Antibody, HRP conjugated

Product Code	CSB-PA004558EB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P41180
Immunogen	Recombinant Human Extracellular calcium-sensing receptor protein (20-612AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Extracellular calcium-sensing receptor (CaR) (CaSR) (hCasR) (Parathyroid cell calcium-sensing receptor 1) (PCaR1), CASR, GPRC2A PCAR1
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	CASR
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.