



HSP90B1 Antibody, FITC conjugated

Product Code	CSB-PA10887C0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P14625
Immunogen	Recombinant Human Endoplasmic protein (27-799AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Endoplasmic (94 kDa glucose-regulated protein) (GRP-94) (Heat shock protein 90 kDa beta member 1) (Tumor rejection antigen 1) (gp96 homolog), HSP90B1, GRP94 TRA1
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	HSP90B1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.