



# IL-33 Antibody, Biotin conjugated

Product Code	CSB-PA011656DD01PI
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	M5B263
Immunogen	Recombinant Pig Interleukin-33 protein (123-276AA)
Raised In	Rabbit
Species Reactivity	Pig
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	IL-33
Immunogen Species	Pig
Research Area	Others
Target Names	IL-33