




BCHE Antibody

Product Code	CSB-PA002604LA01DO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P32750
Immunogen	Recombinant Dog Cholinesterase protein (1-141AA)
Raised In	Rabbit
Species Reactivity	Dog
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cholinesterase (EC 3.1.1.8) (Acylcholine acylhydrolase) (Butyrylcholine esterase) (Choline esterase II) (Pseudocholinesterase) (Fragment), BCHE
Immunogen Species	Dog
Research Area	Others
Target Names	BCHE

Image



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA002604LA01DO at dilution of 1:100

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.