



SFTPb Antibody, Biotin conjugated

Product Code	CSB-PA18055D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P15782
Immunogen	Recombinant Pig Pulmonary surfactant-associated protein B protein (1-79AA)
Raised In	Rabbit
Species Reactivity	Pig
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Pulmonary surfactant-associated protein B (SP-B) (8 kDa protein) (Pulmonary surfactant-associated proteolipid SPL(Phe)), SFTPb, SFTP3
Immunogen Species	Pig
Research Area	Others
Target Names	SFTPb
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.